

Safer Phenol Testing

An important enhancement to Phenol CHEMets & VACU-vials from CHEMetrics

Cleaner. Safer. Same trusted results.

What's New

S-8000 Phenol Activator Solution

Potassium ferricyanide is now provided as a liquid accessory solution.

- ✓ No more loose crystals on the ampoule tip
- ✓ No change to reagent chemistry
- ✓ No change to shelf life
- ✓ No impact on accuracy or performance



Included and Accessories

- Each kit includes 1 bottle of **S-8000 Activator**
- Refill available: **Accessory Pack A-8000**
- New SDS available: **V8000SS**



V-3000 Software Update

Programs **152** and **153** have been updated to reflect new kit code numbers.

- No changes to calibration curve
- No changes to analytical performance
- Test results remain unchanged
- Administrative update only
- No impact on instrument operation

Discontinued Phenol Kits				New Phenol Kits				
Kit Code	Refill	Visual / Instrument	Ferricyanide delivery	Kit Code	Refill	Visual / Instrument	Ferricyanide delivery	Dosing amount
K-8012	R-8012	Visual	Stir to completely dissolve tip-coating of crystals in sample	K-8010	R-8010	Visual	S-8000 Phenol Activator Solution dosed into the sample by drops	5 drops
K-8020D	R-8012	Visual		K-8010D	R-8010	Visual		5 drops
K-8003	-	Instrument		K-8043		Instrument		10 drops
K-8023	-	Instrument		K-8053		Instrument		10 drops

Why it Matters

- Phenol, even at low levels, can indicate contamination or cause taste and odor issues in drinking and environmental waters.
- This update improves phenol testing by:
- Reducing chemical exposure
- Eliminating loose crystals
- Increasing consistency in real-world use