

# Reliable Tools for Crop Testing Labs

BrandTech® Scientific offers the lab tools you need to simplify your pesticide, fungal and potency testing

## The Standard in Bottletop Dispensers

### Dispensette® S dispensers



- Rugged: no seals to wear out
- Suitable for with most reagents
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 121 °C
- Adapters for common bottle threads included
- For trace analysis in ppb, ppt range choose the Dispensette® S Trace Analysis.



Dispensette® S	Without recirculation valve		With recirculation valve	
	Cat. No.	2021 List Price	Cat. No.	2021 List Price
<b>Dispensette® S Digital</b>				
0.1-1	4600310	\$534.80	4600311	\$568.00
0.2-2	4600320	534.80	4600321	568.00
0.5-5	4600330	534.80	4600331	568.00
1-10	4600340	534.80	4600341	568.00
2.5-25	4600350	722.20	4600351	755.40
5-50	4600360	732.80	4600361	766.00
<b>Dispensette® S, Analog-adjustable</b>				
0.1-1	4600100	468.80	4600101	490.00
0.2-2	4600120	468.80	4600121	490.00
0.5-5	4600130	468.80	4600131	490.00
1-10	4600140	468.80	4600141	490.00
2.5-25	4600150	656.20	4600151	677.60
5-50	4600160	672.20	4600161	694.40
10-100	4600170	1,047.20	4600171	1,069.60

*Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.*

Volume, mL	Valve Spring	Without recirculation valve		With recirculation valve	
		Cat. No.	2021 List Price	Cat. No.	2021 List Price
<b>Dispensette® S Trace Analysis, Analog-adjustable</b>					
1-10	Platinum-iridium	4640040	\$1,168.60	4640041	\$1,201.60
1-10	Tantalum	4640240	1,168.60	4640241	1,201.60

## Plant resin extraction process is simplified with VACUUBRAND® vacuum pumps

### PC3001 VARIO® select

**NEW** VARIO® vacuum with the VACUU-SELECT® controller. Automatic control of all common lab applications with touch-screen control!

- Intuitive control
- Touch screen navigation
- Preloaded lab applications
- Fully automatic evaporations
- Drag and drop editing
- Chemical resistant

Description	Cat. No.	2021 List Price
PC 3001 VARIO® select, US plug	20700203	\$7,995.00



**Lab Rats Trust  
BrandTech!**

**BRANDTECH®**  
SCIENTIFIC, INC.

Toll Free: 888-522-2726  
[www.brandtech.com](http://www.brandtech.com)

# Reliable Tools for Crop Testing Labs

## Comfortable, Convenient Meniscus Control

### macro™ pipette controller



- Replaces old-fashioned pipette bulbs with simple, one-finger control
- Rapid or sensitive filling; fast or drop-wise dispensing
- Fits most pipettes 0.1-200 mL




Description	Cat. No.	2021 List Price
<b>macro™ Pipette Controller</b>		
Gray	26200	\$88.80
Blue	26202	88.80
Magenta	26203	88.80
Green	26201	88.80
<b>accu-jet® pro Pipette Controller</b>		
Dark Blue	26330	426.20
Magenta	26331	426.20
Green	26332	426.20
Royal Blue	26333	426.20

## BRAND® Real Time PCR (qPCR)

### Designed to optimize results for qPCR

- DNase-, DNA- and RNase-free
- Suitable for use in common thermal cyclers
- Extra-thin wall dimensions provide optimal thermal transfer and short cycle times



Description	Qty./Pack	Cat. No.	2021 List Price
<b>White PCR Plates (all suitable for qPCR)</b>			
96-well plate, non-skirted, low-profile	50	781367	\$235.60
96-well plate, non-skirted, standard profile	50	781369	235.60
96-well plate, semi-skirted, standard profile	50	781376	235.60
96-well plate, semi-skirted, low-profile	50	781372	235.60
96-well plate, semi-skirted, low-profile, raised skirt	50	781374	235.60
96-well plate, skirted, low-profile	50	781378	381.80
96-well plate, low-profile, for Roche LightCycler®	50	781364	314.00
96-well plate, low-profile, for Roche LightCycler®, with fifty sheets of qPCR sealing film	50	781365	414.00
384-well plate, full skirt, for Roche LightCycler®	50	781358	595.60
<b>qPCR 8-Strip Caps &amp; Tubes</b>			
8 x 0.2 mL tube strip, white	125	781325	170.00
8-strip flat caps, Clear (qPCR)	125	781334	46.00

## Lightweight and Rugged

### Transferpette® S Adjustable Volume Single Channel Pipettes

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Slim design reaches bottom of 15 mL conicals.
- Comfortably fits most hands.
- Accepts tips from most quality manufacturers
- Fixed volume single channels also available



Capacity, µL	Cat. No.	2021 List Price
0.1-1	705868	\$297.00
0.1-2.5	705869	297.00
0.5-10	705870	297.00
2-20	705872	297.00
5-50	705873	297.00
10-100	705874	297.00
20-200	705878	297.00
100-1000	705880	297.00
500-5000 (0.5-5 mL)	705882	297.00
1000-10000 (1-10 mL)	705884	297.00

## BRAND® Class A, USP, Volumetric Glassware

- Accurate measurement instruments for precision analysis
- USP volumetric instruments comply with Class A error limits required by the U.S. Pharmacopeia
- Individually calibrated and supplied with USP batch certificate
- Durable graduations and long service life



Description	Available volume ranges
BLAUBRAND® bulb pipettes, 1 mark, Class A, USP, Certified	0.5 mL to 100 mL
BLAUBRAND® graduated pipettes, Class A, USP, Certified, cotton plug	1 mL to 10 mL
BLAUBRAND® volumetric flask, Class A, USP, Certified, Boro 3.3, glass stopper	5 mL to 2000 mL
BLAUBRAND® Amber volumetric flask, Class A, USP, Certified, Boro 3.3, glass stopper	5 mL to 1000 mL
BLAUBRAND® graduated cylinder, tall form, Class A, USP, Certified, Boro 3.3	5 mL to 2000 mL



**Lab Rats Trust BrandTech!**

**BRANDTECH®**  
SCIENTIFIC, INC.

Toll Free 888-522-2726  
[www.brandtech.com](http://www.brandtech.com)