


Reliable Tools

for Vape and E-liquid Laboratories

The Standard in Bottletop Dispensers

Dispensette® S dispensers

- Rugged: no seals to wear out!
- Compatible with most reagents, including nicotine
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 250 °F (121 °C).
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing



| Dispensette® S | Without recirculation valve | | With recirculation valve | |
|----------------|-----------------------------|-----------------|--------------------------|-----------------|
| | Cat. No. | 2021 List Price | Cat. No. | 2021 List Price |

Dispensette® S Digital

| | | | | |
|--------|---------|----------|---------|----------|
| 0.1-1 | 4600310 | \$534.80 | 4600311 | \$568.00 |
| 0.2-2 | 4600320 | 534.80 | 4600321 | 568.00 |
| 0.5-5 | 4600330 | 534.80 | 4600331 | 568.00 |
| 1-10 | 4600340 | 534.80 | 4600341 | 568.00 |
| 2.5-25 | 4600350 | 722.20 | 4600351 | 755.40 |
| 5-50 | 4600360 | 732.80 | 4600361 | 766.00 |

Dispensette® S, Analog-adjustable

| | | | | |
|--------|---------|----------|---------|----------|
| 0.1-1 | 4600100 | 468.80 | 4600101 | 490.00 |
| 0.2-2 | 4600120 | 468.80 | 4600121 | 490.00 |
| 0.5-5 | 4600130 | 468.80 | 4600131 | 490.00 |
| 1-10 | 4600140 | 468.80 | 4600141 | 490.00 |
| 2.5-25 | 4600150 | 656.20 | 4600151 | 677.60 |
| 5-50 | 4600160 | 672.20 | 4600161 | 694.40 |
| 10-100 | 4600170 | 1,047.20 | 4600171 | 1,069.60 |

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

Flexible discharge tube

| | Cat. No. | 2021 List Price |
|--|----------|-----------------|
|--|----------|-----------------|

For Nominal Volume, mL. PTFE, coiled, length 800 mm (= 31.5") with safety handle.

For Dispensette® S and Dispensette® S Organic

| | | |
|-------------|--------|----------|
| 2, 5, 10 | 708132 | \$200.40 |
| 25, 50, 100 | 708134 | 215.20 |

For seripettor*

| | | |
|-------|--------|--------|
| 2, 10 | 704522 | 135.00 |
|-------|--------|--------|

*Not suitable for use with peroxides due to platinum-iridium valve spring.

A Rugged, Economical Tool for Routine Dispensing

seripettor®

- Dispenses most non-aggressive liquids such as propylene glycol and vegetable glycerin
- Convenient cleaning and maintenance
- Notched volume selector allows for quick and exact volume adjustments
- Comes with thread adapters to fit most common lab bottles



| Volume, mL | Subdivision, mL | A* < ± % | CV* ≤ % | Cat. No. | 2021 List Price |
|------------|-----------------|----------|---------|----------|-----------------|
| 0.2-2 | 0.04 | 1.2 | 0.2 | 4720120 | \$140.40 |
| 1-10 | 0.2 | 1.2 | 0.2 | 4720140 | 140.40 |
| 2.5-25 | 0.5 | 1.2 | 0.2 | 4720150 | 178.00 |

* The value of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20 °C/68° F) and with smooth operation.

A*=Accuracy, CV*=Coefficient of Variation

Serial dispensing

The optional flexible discharge tube with safety handle speeds serial dispensing tasks, and permits fast and precise dispensing even into narrow test tubes. Functions of recirculation valve and safety discharge tube are fully maintained with the flexible discharge tube.



Lab Rats Trust BrandTech!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726

www.brandtech.com

Reliable Tools for Vape and E-liquid Laboratories

Ethylene-acrylate Coated Glass Amber Bottles

The ethylene-acrylate coating envelopes the glass bottle like a protective skin. Should the bottle break, the splintering hazards are reduced considerably.

| | Cat. No. | 2021 List Price |
|--|----------|-----------------|
| Threaded Bottle, PP screw cap, LDPE pouring ring, 250 mL, 33 mm | 704004 | \$28.20 |
| Threaded Bottle, PP screw cap, LDPE pouring ring, 500 mL, 33 mm | 704006 | 38.20 |
| Threaded Bottle, PP screw cap, LDPE pouring ring, 1000 mL, 45 mm | 704008 | 49.40 |
| Bottle Stand, PP, 225 mm x 160 mm x 325 mm | 704275 | 320.00 |




Class A Graduated Cylinders, PMP, Certified

| Volume mL | Tolerance +/-mL | Grad mm | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|--|-----------------|---------|-----------|-------------|-----------|----------|-----------------|
| Class A Graduated Cylinders, PMP, Certified, red screened overprint | | | | | | | |
| 10 | 0.10 | 0.20 | 145 | 15 | 2 | V64614 | \$24.40 |
| 25 | 0.25 | 0.50 | 170 | 22 | 2 | V64714 | 25.60 |
| 50 | 0.50 | 1.00 | 200 | 27 | 2 | V64814 | 26.80 |
| 100 | 0.50 | 1.00 | 250 | 33 | 2 | V64914 | 36.20 |
| 250 | 1.00 | 2.00 | 315 | 44 | 2 | V65014 | 56.40 |
| 500 | 2.50 | 5.00 | 360 | 58 | 1 | V65114 | 45.60 |
| 1000 | 5.00 | 10.00 | 440 | 69 | 1 | V65214 | 66.00 |
| 2000 | 10.00 | 20.00 | 535 | 97 | 1 | V65414 | 106.60 |

| Volume mL | Tolerance +/-mL | Grad mm | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|--|-----------------|---------|-----------|-------------|-----------|----------|-----------------|
| Class A Graduated Cylinders, PMP, Certified, molded graduations | | | | | | | |
| 10 | 0.10 | 0.20 | 145 | 15 | 2 | V64604 | 19.20 |
| 25 | 0.25 | 0.50 | 170 | 22 | 2 | V64704 | 19.20 |
| 50 | 0.50 | 1.00 | 200 | 27 | 2 | V64804 | 20.20 |
| 100 | 0.50 | 1.00 | 250 | 33 | 2 | V64904 | 27.80 |
| 250 | 1.00 | 2.00 | 315 | 44 | 2 | V65004 | 42.60 |
| 500 | 2.50 | 5.00 | 360 | 58 | 1 | V65104 | 43.60 |
| 1000 | 5.00 | 10.00 | 440 | 69 | 1 | V65204 | 64.00 |
| 2000 | 10.00 | 20.00 | 482 | 97 | 1 | V65304 | 104.40 |

- Glass-clear tall-form cylinders calibrated "to contain" Class A per DIN 12681 and ISO 6706
- Red-screened graduations for legibility
- Imprinted lot and year; supplied certificate includes mean test value and lot
- Autoclaving at 121 ° C will not cause permanent exceeding of tolerance limits



Class A, USP, Certified, Graduated Cylinders, Borosilicate Glass

| Volume mL | Subdivision mL | Error limit +/-mL | Height mm | Qty./Pack | Cat. No. | 2021 List Price |
|--|----------------|-------------------|-----------|-----------|----------|-----------------|
| BLAUBRAND® graduated cylinder, tall form, class A, USP, certified, Boro 3.3 | | | | | | |
| 5 | 0.1 | 0.05 | 115 | 2 | 32805 | \$35.00 |
| 10 | 0.2 | 0.1 | 140 | 2 | 32808 | 35.00 |
| 25 | 0.5 | 0.25 | 170 | 2 | 32820 | 36.20 |
| 50 | 1 | 0.5 | 200 | 2 | 32828 | 41.20 |
| 100 | 1 | 0.5 | 260 | 2 | 32838 | 43.00 |
| 250 | 2 | 1 | 335 | 2 | 32848 | 82.00 |
| 500 | 5 | 2.5 | 365 | 2 | 32854 | 141.40 |
| 1000 | 10 | 5 | 465 | 1 | 32862 | 99.20 |
| 2000 | 20 | 10 | 505 | 1 | 32864 | 184.20 |

For additional BLAUBRAND USP glassware see www.brandtech.com

BrandTech® offers BLAUBRAND® USP volumetric instruments that comply with the Class A error limits required by the United States Pharmacopeia.



Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice.



Toll Free 888-522-2726

www.brandtech.com