

Easy Adjustments
Through Shaft for quick impeller and batch changes

No Downtime
Continuous viscous mixing without fail

Precision
Toggle display to view RPM or Ncm to track changes in viscosity


Efficiency
No air escapes the totally enclosed rugged metal housing

3000 RPM
Dispersion
Pigment integration and high shear mixing

Emulsion
Low shear high pumping pitched blade

Cool Down
Anchor paddle for high viscosities

20 RPM

Canadian Made 

The Kit:

	Universal Overhead Stirrer	Safety Stand	Fastening Clamp	Impeller Shaft	Dispersion Blade	Pitched Blade Propeller	Anchor Paddle
Caframo Part	BDC3030	A110	A120	A722	A164	A531	U044