

Reliable Tools for Industrial Labs

The Standard in Bottletop Dispensers

Dispensette® S dispensers

- Rugged: no seals to wear out!
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 250 °F (121 °C)
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing
- Models to dispense most solvents and acids, including HF (not listed).



Volume, mL	Cat. No.	Without recirculation valve	With recirculation valve	
		2021 List Price	2021 Cat. No.	2021 List Price
Dispensette® S Digital				
0.1-1	4600310	\$534.80	4600311	\$568.00
0.2-2	4600320	534.80	4600321	568.00
0.5-5	4600330	534.80	4600331	568.00
1-10	4600340	534.80	4600341	568.00
2.5-25	4600350	722.20	4600351	755.40
5-50	4600360	732.80	4600361	766.00
Dispensette® S, Analog-adjustable				
0.1-1	4600100	468.80	4600101	490.00
0.2-2	4600120	468.80	4600121	490.00
0.5-5	4600130	468.80	4600131	490.00
1-10	4600140	468.80	4600141	490.00
2.5-25	4600150	656.20	4600151	677.60
5-50	4600160	672.20	4600161	694.40
10-100	4600170	1,047.20	4600171	1,069.60

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

Faster, Easier Titrations with Class A Precision

Titrette® bottletop burette




- Meets Class A standards for glass burettes!
- Popular for sulfur dioxide (SO₂) analyses, acid-base titrations, solvent analyses and more
- Suitable for NaOH, KOH, KI, H₂SO₄, KMnO₄, Iodine and more
- Eliminates meniscus reading errors
- Models available with direct data transmission to PC

Description	Cat. No.	2021 List Price
10 mL (standard)	4760141	\$1,092.00
25 mL (standard)	4760151	1,092.00
50 mL (standard)	4760161	1,092.00

Precision analysis requires measuring instruments with consistently high precision. BLAUBRAND® volumetric instruments represent the highest degree of technical perfection. They are made of glass which is highly resistant to nearly all chemical substances.

- Volumetric instruments are individually calibrated (USP batch certificate supplied)
- High quality inks and annealing process provide durable graduations and long service life
- Thermal stress in the glass blanks is eliminated prior to adjustment so there are no permanent volume changes after heating, up to 250 °C



BLAUBRAND® graduated cylinder, tall form, Class A, USP, Certified, Boro 3.3

Volume (mL)	Height (mm)	Qty./Pack	Cat. No.	2021 List Price
5	115	2	32805	\$35.00
10	140	2	32808	35.00
25	170	2	32820	36.20
50	200	2	32828	41.20
100	260	2	32838	43.00
250	335	2	32848	82.00
500	365	2	32854	141.40
1000	465	1	32862	99.20
2000	505	1	32864	184.20

Convenient & Economical

Disposable PS and PMMA Cuvettes

- Consistent reliable results—manufactured from the highest quality virgin resins
- Accommodate most applications - available in macro and semi-micro sizes
- Easy to use—no more cleaning



Lab Rats Trust BrandTech!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726
www.brandtech.com

Reliable Tools for Industrial Labs

Precise, Slim, Light & Dependable Stepper

HandyStep® S Repeating Pipette



The HandyStep® S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® touch and HandyStep® touch S models also available

Description	Cat. No.	2021 List Price	
HandyStep® S Repeating Pipette	705110		\$402.80
HandyStep® touch Repeating Pipette	705200		1,055.00
PD-Tip™ // size		Non-sterile tips	BIO-CERT® Pure*
0.1 mL (100 tips)	705700	\$125.20	705730 \$171.60
0.5 mL (100 tips)	705702	125.20	705732 171.60
1.0 mL (100 tips)	705704	125.20	705734 171.60
1.25 mL (100 tips)	705706	125.20	705736 171.60
2.5 mL (100 tips)	705708	125.20	705738 171.60
5 mL (100 tips)	705710	125.20	705740 171.60
10.0 mL (100 tips)	705712	125.20	705742 171.60
12.5 mL (100 tips)	705714	125.20	705744 171.60
25 mL (50 tips)	705716	84.80	—
25 mL (25 tips)	—	42.40	705746 52.60
50 mL (25 tips)	705718	42.00	705748 52.60
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	169.60	—

*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

Comfortable, Convenient Meniscus Control

macro™ Pipette Controller



- Replaces old-fashioned pipette bulbs with simple, one-finger control
- Rapid or sensitive filling; fast or drop-wise dispensing
- Fits most pipettes 0.1-200 mL

Description	Cat. No.	2021 List Price
macro™ pipette controller		
Gray	26200	\$88.80
Blue	26202	88.80
Magenta	26203	88.80
Green	26201	88.80

Lightweight and Rugged

Transferpette® S Adjustable Volume Single Channel Pipettes



- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Slim design reaches bottom of 15 mL conicals.
- Comfortably fits most hands.
- Accepts tips from most quality manufacturers
- Multichannel and fixed volume single channel pipettes also available

Capacity, µL	Cat. No.	2021 List Price
0.1-1	705868	\$297.00
0.1-2.5	705869	297.00
0.5-10	705870	297.00
2-20	705872	297.00
5-50	705873	297.00
10-100	705874	297.00
20-200	705878	297.00
100-1000	705880	297.00
500-5000 (0.5-5 mL)	705882	297.00
1000-10000 (1-10 mL)	705884	297.00



Compact Vacuum Pump for Filtration

VACUUBRAND® ME 1 Vacuum Pump



VACUUBRAND® ME 1 is a compact, high-performance, easy-to-use vacuum pump which is ideal for single and multiple filtrations. An optional vent control valve with dial gauge is available for variable adjustment of pumping speed.

- Small footprint conserves valuable bench space
- Whisper quiet performance & low vibration minimize background noise in the lab
- PTFE diaphragms and valves offer good chemical resistance
- High performance 14 lpm pumping speed & 100 mbar ultimate vacuum

Model	Cat. No.	2021 List Price
ME 1 Vacuum Pump	20721003	\$911.00
Vacuum regulation valve with manometer for ME 1	20696842	194.80

Need Lab Vacuum?

Use the [BrandTech® Vacuum Selection Guide](#) to find the best vacuum pump for a wide range of applications! Go to www.brandtech.com.