




Shamrock Tape & Labels for Water Testing

PRESERVATIVE IDENTIFICATION LABELS


Circle Labels

Identify container contents on cap and/or sides with eye-catching 1" dot labels. Permanent adhesive holds securely throughout use. Labels are supplied 1,000 per roll.



Tamper Evident Labels

Place identifying center circle of label on cap and smooth ends down each side of the container. Unique material features serrated edges that shred to indicate tampering if removal is attempted. Labels measure 1 3/8" x 9" and are supplied 100 labels per package.



TAMPER RESISTANT LABEL & TAPE

Label and tape feature destructible material and permanent adhesive. When removal is attempted, material tears, leaving portions of the material in place to immediately signal tampering.

PRE-CUT LABELS

EL-11 1 3/8" x 6 5/8"
100 Individual Labels per Pack
Minimum Order 6 Packs



TAPE LABELS

EL-110-4 1 3/8" x 108"
Minimum Order 10 Rolls



COLLECTION LABEL


This label's specially treated material will hold securely, even after contact with water, yet can be easily removed for container reuse. Labels are supplied 500 per roll. Call for information on labels imprinted with your facility name.

Collection Date _____
Collection Time _____
Collector _____
Lab Sample # _____
Address _____
1st Draw _____ 2nd Draw _____
Preservative _____

WT-50 2 1/4" x 4"

BIOHAZARD WARNING LABELS

Labels feature bright, red fluorescent color and permanent, pressure-sensitive adhesive to hold securely throughout use.



SBH-101 8" x 10" Pressure-sensitive Package of 10



SBH-5 3 1/2" x 2 1/4" Package of 100



SBH-1 3" x 2 1/4" Roll of 500

STANDARD BLANK LABELS

Smudgeproof white labels offered in your choice of standard removable adhesive for use on clean, smooth surfaces, or in premium Super-Stick removable material which is perfect for small, irregular and hard to label surfaces. Supplied 1,000 labels per roll.




Size	Standard	Super-Stick
1" x 3"	ADS-1	ADSF-1
1" x 3 1/2"	ADS-2	ADSF-2
1 1/4" x 3"	ADS-3	ADSF-3
1 1/4" x 3 1/2"	ADS-4	ADSF-4


Size	Standard	Super-Stick
1 1/2" x 3"	ADS-5	ADSF-5
1 1/2" x 3 1/2"	ADS-6	ADSF-6
2" x 3"	ADS-7	ADSF-7
2 1/2" x 4"	ADS-8	ADSF-8

BLANK THERMAL LABELS

A few examples are illustrated below. Contact Shamrock for a complete listing of thermal label sizes and materials available.




THERM-563R Thermal Direct
THERM-563R-T Thermal Transfer
1 3/4" x 3" 3" Core 1,800/roll



THERM-358R Thermal Direct
THERM-358R-T Thermal Transfer
2" x 4" 3" Core 1,600/roll



THERM-1791R-DO
1 1/4" x 2 1/4"
1,000 per roll



THERM-1792R-DO
4" x 2 1/4"
250 per roll

CONSECUTIVE NUMBER LABELS

Labels feature glove friendly adhesives that stick securely to plastic, glass and metal containers. Contact Shamrock for complete information on the variety of label and number options available, and for ordering assistance.

Tape Numbers


Available in 19 colors, tape labels resist moisture and temperature extremes yet are easily removed, may be transferred to another surface. Labels measure 1/2" x 1 1/4".



SCR-10 (500 labels per roll) SCR-20 (2000 labels per roll)

Pre-Cut Label Numbers


Choose tablet or roll form. Permanent adhesive holds fast on hard to label surfaces. Labels measure 1/2" x 1 1/4" and are available in 19 colors.



SCS-12 (Tablet) SCR-14 (Roll)

Bar Coded Numbers

Available in tablets or rolls and in a variety of standard formats, white labels feature sharply printed bar codes and eye-readable numbers.




Tablet Form

MICROSCOPE SLIDE LABELS


Standard Labels

Supplied 1,000 labels per package or roll. Labels measure 1 5/16" x 1 5/16".



Thermal Labels

Available for thermal direct or thermal transfer printing.



THERM-1213R Thermal Direct
THERM-1213R-T Thermal Transfer
THERM-1213OPR Thermal Transfer Opaque Polyester
1 5/16" x 1 5/16" 3" Core 12,800/roll



THERM-1839R-D0
Thermal Direct
1 5/16" x 1 5/16" 1" Core 750/roll

Offered in a variety of formats, all feature permanent adhesive to hold fast during use and storage.

Pinfed Labels

Supplied 30,000 labels per carton.



BDP-1 1 5/16" x 1 5/16" Fan folded

Mylar Overlay Label


Protect handwritten or printer generated information. Labels are supplied 500 per roll.



ME-24 1" x 1"

Laser Labels

Supplied on 8 1/2" x 11" sheets with 5,000 labels per package.



SLS-LASER-RC 1 5/16" x 1 5/16"
80 labels per sheet



Shamrock

SHAMROCK SCIENTIFIC SPECIALTY SYSTEMS, INC.
34 DAVIS DR., BELLWOOD, IL 60104

Water Testing Labels




VISIT SHAMROCK AT
www.ShamrockLabels.com



PRSR STD
U.S. POSTAGE
PAID
Carol Stream, IL
Permit No. 925

Fold


CALL 1-800-323-0249 • FAX 1-800-248-1907 • EMAIL sales@ShamrockLabels.com • SHOP www.ShamrockLabels.com 913



Size	Catalog No.
1/4"	SDL-25
7/16"	SDL-43
1/2"	SDL-50
3/4"	SDL-75
1"	SDL-100
1 1/4"	SDL-125

- 6 Sizes • 31 Colors
- Perfect for Color Coding by test, collector, draw
- Your choice of permanent or removable adhesive
- Available in sheet form or on rolls in dispenser boxes
- May be imprinted to your specifications

DOT LABELS



Width	Cat. No.	500"	2160"
1/2"	ST-12	ST-120	ST-120
3/4"	ST-34	ST-340	ST-340
1"	ST-10	ST-100	ST-100
1 1/2"	ST-11	ST-110	ST-110
2"	ST-20	ST-200	ST-200

- 5 Widths • 2 Roll Lengths • 31 Colors
- Resists Moisture, Solvents, Temperature Extremes
- Easily Removed from Glass, Plastic, Metal or Wood
- Dispensers available for single or multiple rolls
- May be imprinted to your specifications

LABELING TAPE

Shamrock Labeling Tape & Dot Labels