


Reliable Tools for Biopharm Lab



BRANDplates® microplates

For Cell Culture, Immunoassay & General Use

- 9 different surface treatments
- Transparent, white, black, colored with clear bottom
- C, F, V and U-shaped wells
- Deep well plates
- Insert System for advanced cell culture
- Barcoding on request
- Samples available



BRAND® offers a comprehensive line of microplates for general, cell culture and immunology applications. The line includes plates with distinctive physical treatments that eliminate refrigeration and shelf-life concerns. All are manufactured according to ANSI-SLAS 1 to 4-2004. Visit www.brandtech.com for additional information.

Dispensette® S dispensers

The Standard in Bottletop Dispensers

- The 45 mm standard thread, plus included adapters, fit most common lab bottles
- The instrument is easy to disassemble for cleaning
- Valves are replaceable for simple, economical service
- New valve design eliminates seals for trouble-free operation
- Dispensette® S is autoclavable at 121 °C
- Easy to calibrate and adjust in order to comply with ISO 9001 and GLP guidelines



Dispensette® S	Without recirculation valve		With recirculation valve	
	Cat. No.	2021 List Price	Cat. No.	2021 List Price
Dispensette® S Digital				
0.1-1	4600310	\$534.80	4600311	\$568.00
0.2-2	4600320	534.80	4600321	568.00
0.5-5	4600330	534.80	4600331	568.00
1-10	4600340	534.80	4600341	568.00
2.5-25	4600350	722.20	4600351	755.40
5-50	4600360	732.80	4600361	766.00
Dispensette® S, Analog-adjustable				
0.1-1	4600100	468.80	4600101	490.00
0.2-2	4600120	468.80	4600121	490.00
0.5-5	4600130	468.80	4600131	490.00
1-10	4600140	468.80	4600141	490.00
2.5-25	4600150	656.20	4600151	677.60
5-50	4600160	672.20	4600161	694.40
10-100	4600170	1,047.20	4600171	1,069.60

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

The BVC family of liquid aspiration systems

The BVC family of fluid aspiration systems deliver compact and quiet solutions for cell culture work. They provide the flexibility to accommodate your decontamination protocol whether you autoclave or disinfect in the bottle.

- Reliable: VACUUBRAND® pump* for long service life
- Power & Control: Suitable for two users, with adjustable* vacuum level
- Robust: "G" models designed for in-bottle bleach decontamination
- Safe: aerosol protection; models with contact-free level sensor available
- Single and multichannel hand controllers for use with cell culture vessels and multi-well plates

* Control and professional models



Lab Rats Trust BrandTech®!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726 | www.brandtech.com

Reliable Tools for Biopharm Labs

Precise, Slim, Light & Dependable Stepper

HandyStep® S Repeating Pipette



The HandyStep® S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® touch and HandyStep® touch S models also available

Description	Cat. No.	2021 List Price	
HandyStep® S Repeating Pipette	705110		\$402.80
HandyStep® touch Repeating Pipette	705200		1,055.00
PD-Tip™ // size		Non-sterile tips	BIO-CERT® Pure*
0.1 mL (100 tips)	705700	\$125.20	705730 \$171.60
0.5 mL (100 tips)	705702	125.20	705732 171.60
1.0 mL (100 tips)	705704	125.20	705734 171.60
1.25 mL (100 tips)	705706	125.20	705736 171.60
2.5 mL (100 tips)	705708	125.20	705738 171.60
5 mL (100 tips)	705710	125.20	705740 171.60
10.0 mL (100 tips)	705712	125.20	705742 171.60
12.5 mL (100 tips)	705714	125.20	705744 171.60
25 mL (50 tips)	705716	84.80	— —
25 mL (25 tips)	—	—	705746 52.60
50 mL (25 tips)	705718	42.40	705748 52.60
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	169.60	— —

*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

Comfortable, Convenient Meniscus Control

macro™ Pipette Controller



- Replaces old-fashioned pipette bulbs with simple, one-finger control.
- Rapid or sensitive filling; fast or drop-wise dispensing.
- Fits most pipettes 0.1 to 200 mL.

Description	Cat. No.	2021 List Price
macro™ Pipette Controller		
Gray	26200	\$88.80
Blue	26202	88.80
Magenta	26203	88.80
Green	26201	88.80
accu-jet® pro Pipette Controller		
Dark Blue	26330	426.20
Magenta	26331	426.20
Green	26332	426.20
Royal Blue	26333	426.20

accu-jet® pro Pipette Controller



- Provides smooth, intuitive meniscus control.
- Simple, sensitive soft blow-out.
- Contoured grip for comfort.
- Mode switch permits either powered or gravity dispensing.

Lightweight, Robust & Fully Autoclavable

Transferpette® S Adjustable Volume Single Channel Pipettes



Capacity, µL	Cat. No.	2021 List Price
0.1-1	705868	\$297.00
0.1-2.5	705869	297.00
0.5-10	705870	297.00
2-20	705872	297.00
5-50	705873	297.00
10-100	705874	297.00
20-200	705878	297.00
100-1000	705880	297.00
500-5000 (0.5-5 mL)	705882	297.00
1000-10000 (1-10 mL)	705884	297.00

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain.
- Slim design reaches bottom of 15 mL conicals.
- Comfortably fits most hands.
- Accepts tips from most quality manufacturers
- Fixed volume single channels also available.

Transferpette® S -8/-12 Multichannel Pipettes



Capacity, µL	8 Channel	2021 List Price	12 Channel	2021 List Price
0.5-10	705900	\$760.00	705920	\$940.00
5-50	705906	760.00	705926	940.00
10-100	705908	760.00	705928	940.00
20-200	705910	760.00	705930	940.00
30-300	705912	760.00	705932	940.00

Transferpette® S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, reagent reservoir, pack of replacement O-rings. One-year warranty.



Sample Cups for Clinical Analyzers

BRAND® sample cups are made from virgin polystyrene and are designed to fit many common analyzers. Caps and lids are made of polyethylene for a snug fit. Cups and caps meet CE-IVD Directive 98/79



Volume, mL	Capacity, mL	Cat. No.	2021 List Price
Sample Cups for Technicon® Analyzers	1.5	115015	\$550.20
Sample Cups for Technicon® Analyzers	2	115016	590.00
Sample Cups for Technicon® Analyzers	4	115017	415.60
Push Caps (for 115015 or 115016)		115020	58.40
Coulter Counter® with Lids	20	722055	181.00

Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice.

Visit www.brandtech.com for additional information.