

The leaves are falling...
and so is the price on the BeadBlaster™ 24!



Offer Valid Through: 12.31.18
Promo Code: FP2400

BEADBLASTER™ 24

MICROTUBE HOMOGENIZER



\$6495*
Valid thru 12/31/18

~~\$7995~~

- + Extremely powerful mixing action, for lysis, grinding or homogenization
- + Sealed 2.0ml tubes - no cross-contamination
- + Faster & more effective than tissue grinders
- + Ideal for extraction of DNA, RNA and Proteins from all sample types



Zirconium
Triple-Pure™: Acid washed, heat treated & tested free of nuclease

**BUY TWO PACKS OF BEADS,
GET ONE FREE
WITH THE PURCHASE OF A
BEADBLASTER™ 24**

Triple-Pure™ (Molecular Biology Grade) Zirconium Beads

Ordering Information:

Item No.	Description:	Price:
D2400	BeadBlaster™ 24 Microtube Homogenizer* with 24 x 2.0ml tube holder	\$6495
D2400-R5	Optional tube holder, 12 x 5ml	\$345.00
Prefilled Tube Kits:	Include 2.0ml non-skirted screw cap tubes with caps, sealing rings and beads as specified below:	
D1031-01	Standard glass beads, Ø 0.1mm, 50pk.	\$104.50
D1031-05	Standard glass beads, Ø 0.5mm, 50pk.	\$104.50
D1031-10	Standard glass beads, Ø 1.0mm, 50pk.	\$104.50
D1032-01	Triple-Pure, high impact zirconium beads, Ø 0.1mm, 50pk.	\$146.00
D1032-05	Triple-Pure, high impact zirconium beads, Ø 0.5mm, 50pk.	\$146.00
D1032-10	Triple-Pure, high impact zirconium beads, Ø 1.0mm, 50pk.	\$146.00
D1032-15	Triple-Pure, high impact zirconium beads, Ø 1.5mm, 50pk.	\$146.00
D1032-30	Triple-Pure, high impact zirconium beads, Ø 3.0mm, 50pk.	\$146.00
D1032-SK	Starter Kit, 10 ea. of Ø 0.1, 0.5, 1.0, 1.5 & 3.0 Triple-Pure, high impact zirconium	\$146.00
D1032-60	One ceria zirconium oxide satellite, Ø 6.0mm, 50pk.	\$146.00
D1033-28	Stainless steel beads, Ø 2.8mm, 50pk.	\$146.00
D1033-30G	Garnet shards with 1 zirconium bead Ø 6.0mm, 50pk.	\$146.00
D1034-MX	Mixed format, 0.1mm silica, 1.4mm zirconium, 4mm silica, 50pk	\$146.00
D1031-T20	Empty tubes (2.0ml) pack of 1000	\$154.00
D1031-RF	Reinforced empty tubes (2.0ml) pack of 500	\$189.00
D1034-28	5ml tubes (snap cap) with stainless beads, Ø 2.8mm, 50pk	\$195.00

Additional bead types are also available (in both bulk and prefilled kits). Please contact Benchmark for a complete listing.



Scan the QR code to watch the BeadBlaster™ 24 in action!

Technical Data:

Capacity:	24 x 2.0ml tubes
Speed:	4.00 to 7.00 m/s, 0.05 increments
Cycle time:	1 - 90 sec, 1 sec increments
Pause between cycles:	0-2 min, 1 sec increments
Max. cycles per program:	10
Program memory:	up to 50
Accel. / Decel. time:	<2 sec / <2 sec
Noise level:	<68db
Dimensions (WxDxH):	11 x 14.2 x 15.2 in / 28 x36 x38.5 cm
Weight:	55 lb / 25kg
Electrical:	100V to 240V, 50-60Hz, 600W
Warranty:	2 years

HERMLE Z326 Universal Centrifuge



- Compact and versatile
- 17 rotor options with a capacity for 4x100ml in a swing out rotor
- Maximum 18,000rpm/23,524xg
- Refrigerated and non-refrigerated options

Z326-BND	Z326 Universal Centrifuge Bundle includes centrifuge, 4x100ml swing out rotor, 3x15ml buckets (pk of 2) and 2x50ml buckets (pk of 2)	\$6,749.00
Z326-K-BND	Z326 K Universal Refrigerated Centrifuge Bundle (as above)	\$8,849.00
Z326-02MP	2 x 3 microplate rotor	\$2,162.19
Z326-1215	12 x 15ml rotor	\$1,031.87
Z326-3015	30 x 15 ml rotor	\$1,595.05
Z326-0685	6 x 85ml high speed (38°)	\$2,295.00
Z326-0650LS	6 x 50ml conical (6,000rpm)	\$1,135.55
Z326-2420	24 x 1.5/2.0ml rotor	\$611.77
Z326-4420	44 x 1.5/2.0ml rotor	\$1,031.87
Z326-3202	4 x PCR strip rotor	\$1,195.00

Additional rotors, buckets and inserts available and sold separately. Visit www.hermleusa.com for details.

BenchMixer™



- Powerful motor for "instant vortexing"
- Continuous or touch operation
- Q-drive counter balance drive system
- Wide range of accessories

BV1000	BenchMixer™ Vortex Mixer w/cup head	\$239.00
BV1000-COMBO	Combo head for plates and microtubes	\$59.00
BV1000-H500	Horizontal head for 2x50ml tubes	\$85.00

Additional heads available

BEADBUG™ 6

MICROTUBE HOMOGENIZER



- Fast, efficient bead homogenization
- Eliminates cross contamination
- High speed, 3D motion
- Accommodates 2ml or 5ml tubes
- Large variety of beads available

D1036	BeadBug™ 6 Microtube Homogenizer	\$1995.00
D1032-01	Prefilled Tubes, 0.1mm zirconium beads, pk/50	\$146.00
D1032-30	Prefilled Tubes, 3.0mm zirconium beads, pk/50	\$146.00

Additional bead sizes and materials available

INCU-SHAKER™

Mini



- Touch screen control panel
- Compact design - just 11" wide
- Exceptional temperature uniformity
- 19mm orbit
- Optional MAGic Clamp™ accessories for tool free exchange

H10001-M	IncuShaker™ Mini w/flat platform, nonslip mat	\$2895.00
H1000-MR	MAGic Clamp™ Platform for use w/MAGic Clamps	\$209.00
H1000-MR-CMB	Adjustable Flask Clamp for 125, 250 & 500ml flasks	\$37.50
H1000-MR-1550	MAGic Clamp Tube Rack for 30x15ml and 20x50ml	\$129.00

A variety of MAGic Clamp accessories available

MY GEL™

mini



- Complete system: includes gel box, 2 sizes of gel trays, casting stand and integrated power supply
- Quick and easy casting, no tape or spacers
- Ideal for personal use

E1101	myGel™ Mini Electrophoresis System with integral power supply	\$495.00
E1101-SK	myGel Mini Electrophoresis System Starter Kit includes myGel Mini, 100g agarose and mini balance	\$545.00

Digital Hotplate Stirrers



- Large, backlit LCD
- 6.5"x6.5" ceramic work surface
- Control actual sample temperature with optional probe
- 3 models - heat/stir, heat, stir

H3760-HS	Digital Hotplate Magnetic Stirrer	\$349.00
H3760-H	Digital Hotplate	\$229.00
H3760-S	Digital Magnetic Stirrer	\$229.00
H3760-TP	Optional Temperature Probe	\$99.00
H3760-CS	Clamp set for H3760 series, includes rods and clamps	\$68.00

Effective 9/15/2018. Pricing, availability and specifications subject to change without notice



Fox Scientific, Inc.
 8221 East FM 917 • Alvarado, TX 76009
 Phone: (800)369-5524 • Fax: (817)783-3571