

Battery Powered Magnetic Stirrer

Stirring Power – Anywhere, Anytime



Ideal for field use, incubators or anywhere electricity is not convenient or not available. Provides up to 16 days of continuous operation with 2 "D" batteries.

- Mix up to 1 liter
- Speed control knob up to 900RPM
- Use with magnetic stirring bars 1" or smaller
- Polypropylene construction offers durability, chemical resistance and ease of cleaning

Catalog No.	Dimensions	List Price/Each
H37017-0000	9¼ x 5 x 3½"	\$216.15

Labmat™ Rolls & Sheets

White or Biohazard Safety Orange

Working with chemicals, dyes, corrosive chemicals, radioactive solutions, or other liquids? Disposable Labmat™ Bench Liner speeds clean-up and safeguards work surfaces from spills, wear and stains.

- Absorbent upper surface traps up to 420ml/m2 of liquid
- Non-skid polyethylene backing prevents soak through of most solutions
- Easily cut with scissors



Catalog No.	Description	List Price/Each
F24675-0000	White 20" x 50 ft. roll	\$52.00
F24675-1000	Safety Orange for biohazard designated use. 20" x 50ft. roll	63.80
F24675-0020	100 Cut Sheets, (White), 20 x 24"	184.60

Disposal Cartons for Glass

Dispose of Glass Safely

Sturdy corrugated receptacle comes complete with 2mil polyethylene bag to contain glass fragments and liquid spills. Snug fitting lid has a flap opening for inserting waste. When full, pull the flap closed and the entire carton is ready to dispose. Labeled in English, French, German and Spanish, they are packaged flat to save space and shipping costs and are quickly assembled.



Catalog No.	Model	Dimensions (L x W x H)	Shipping Weight	List Price/Pack
F24653-0001	Floor	12 x 12 x 27"	6.8kg (15lb)	\$73.65/6
F24653-0002	Benchtop	8 x 8 x 10"	3.2kg (7lb)	52.80/6

Spinpak® Magnetic Stirring Bar Assortment with Restrainer

Large to Small Volumes, the Spinpak® Assortment has You Covered.

10-pack of the most popular sizes and colors of Spinbar® octagonal magnetic stirring bars.

- Octagonal shape provides greater surface area for added turbulence compared to smooth cylindrical bars
- Integral pivot ring aids in reducing friction and chattering
- Spinbar restrainer fits in one hand and can be used to hold bars up to 3" long against vessel walls, eliminating splashing hazards when pouring



Catalog No.	Spinpak® Sizes	List Price/Each
H37158-0010	2 each – ½ x ½", ½ x ¾", 1 x ¾", 1½ x ¾", 2 x ¾"	\$63.75

Safety Waste Jug

Collect & Transport Spent Solvents

This high-density polyethylene container is ideal for temporary storage of spent solvents and other liquid waste. Screw-on funnel top makes it easy to safely collect materials and a snap-on cap seals the funnel top to retain fumes. Funnel top can be replaced with an included screw-on top for safe transit.



Catalog No.	Description	Dimensions	List Price/Each
H11919-0000	22 liters (5.5 gallons)	9¾ x 10¼ x 14¼"	\$126.30

Acid/Solvent Bottle Carriers

Polyethylene carriers feature wide ribs in the sides and bottom to support and cushion bottles and their contents. Choose from two styles of openings: Cone style has a 5/8" opening for the insertion of tubing so liquids can be transferred to the container without pouring. The truncated style permits pouring without removal of the lid. Designed to hold standard reagent bottles.



Catalog No.	Style	Holds Bottle	List Price/Each
F16957-0000	Cone	2.5 liter (5 pints)	\$71.55
F16958-0000	Cone	4 liter (1 gallon)	84.95
F16959-0000	Truncated	2.5 liter (5 pints)	71.55
F16960-0000	Truncated	4 liter (1 gallon)	84.40

Lit. No. 146 - 9/11



SAMPLE HANDLING & BEYOND

Gather it, organize it, store it, dispense it, process it, dispose it . . .
Do it faster, easier & safer with Scienceware® tools and supplies

Prevent Accidental Misuse, Keep Liquids in Clearly Marked Bottles!

Right-to-Know, Vented 4-Color Wash Bottles



Preprinted safety-vented wash bottles for commonly used chemicals in the workplace feature a patented valve that prevents pressure build-up of volatile solvents and eliminates spontaneous discharge. Vented bottles can be dispensed from an upright or inverted position. U.S. Patent No. 5,125,543. 8oz and 32oz sizes and non-vented bottles also available for many liquids. Pre-printed Bottles feature:

- Chemical Name & Formula
- Full Color NFPA Diamond
- Health Hazards
- CAS#
- Suggested Protective Clothing and Equipment

Catalog No.	Label	Size	List Price/Pack
F11816-0050	Assortment †	500ml (16oz)	57.05/6
F11816-0001	Acetone	500ml (16oz)	38.10/4
F11816-0004	Distilled Water	500ml (16oz)	38.10/4
F11816-0019	Ethanol	500ml (16oz)	38.10/4
F11816-0007	Ethyl Acetate	500ml (16oz)	38.10/4

*Note: Use LYOB only with solvents compatible with LDPE. **Sodium Hypochlorite bottle is white LDPE for light-protection between 230 and 450nm.
†Assortment is 6pcs., 4 assortments per case and includes one each of Acetone, Bleach, Distilled Water, Ethanol, Isopropanol and LYOB.

Bottle Carrier

Portable Convenience

Made of ¼" thick high-density polyethylene, this carrier is 14 x 7 x 6" high with a 1" diameter handle. Use with or without the divider (sold separately) which creates 6 compartments measuring 3¾" x 3¾" and will hold 8 oz and 16oz bottles.

Catalog No.	Description	List Price/Each
F16963-0000	Without partitions	\$171.40
F16962-0006	6 compartment divider insert	50.95

Wash Bottle Carrier

Easily Transport 8 or 16oz. Bottles

Lightweight, polypropylene carrier holds six, 8 or 16 oz. wash bottles in 2.9" diameter openings. Autoclavable.

Catalog No.	Description	List Price/Each
F11640-0001	Wash Bottle Carrier, 11¼ x 7¾ x 8¾"	\$37.75



Safety-Vented Patented Valve

Wide Mouth (1½") for Easy Filling



Back View

Safety Vent

- No Spontaneous Discharge
- Dispenses Upright or Inverted





Precisionware®
Low-Density
Polyethylene Bottles
 FDA Compliant Material;
 Leakproof Design

- Choose from narrow- or wide-mouth design
- Perfect for sampling and general storage of liquids, powders or viscous fluids
- Translucent, flexible LDPE with PP closures

Narrow-Mouth Bottles

Catalog No.	Capacity	Closure	List Price/Pack
F10621-0015	125ml (4 oz)	28mm	\$21.45/12
F10621-0016	250ml (8oz)	28mm	24.40/12
F10621-0007	500ml (16oz)	28mm	32.45/12
F10621-0008	1000ml (32oz)	38mm	29.75/6

Wide-Mouth Bottles

Catalog No.	Capacity	Closure	O.D.	List Price/Pack
F10626-0005	125ml (4oz)	38mm	1.9"	\$23.65/12
F10626-0006	250ml (8oz)	45mm	2.3"	31.65/12
F10626-0007	500ml (16oz)	53mm	2.84"	44.50/12
F10626-0018	1000ml (32oz)	53mm	3.56"	41.00/6





Where You Want It – When You Need It

PLUM Instant Access Eye Wash Stations

Easy-to-Mount work station units reduce critical time and distance to first aid, making them a perfect supplement to fixed stations. Sterile 0.9% Saline effectively rinses debris and other foreign matter from the eyes and PLUM exclusive pH Neutralizing Wash is ideal for accidents involving acids and alkalis, effectively neutralizing either in 20 seconds or less.

- Built-in flexible eye cups; quick twist open
- Bottles can be used upright or inverted
- Easy-to-follow pictograms instructions
- Stations easily mount to walls and other flat surfaces; hardware is included
- Great for laboratories, manufacturing, classrooms
- Covered stations are ideal for high traffic/dusty areas



Accidents can and do happen to both eyes at the same time. PLUM Duo bottles give you the ability to treat both eyes simultaneously.



Rinse times: pH Neutralizing Wash – approximately 2 min/bottle, Saline Eye Wash – approximately 5 min./bottle.

Catalog No.	Description	List Price/Each
A. F24880-4003	Open Station - 200ml (6.8oz) pH Neutralizing and 500ml (16oz) 0.9% Saline Eye Wash, 9 x 12"	\$57.20
B. F24880-5003	Covered Station - 200ml (6.8oz) pH Neutralizing and 500ml (16oz) 0.9% Saline Eye Wash, 8 3/4 x 11"	67.60
C. F24880-8403	DUO Open Station - 500ml (16oz) pH Neutralizing and 1000ml (32oz) 0.9% Saline, 12 x 15"	\$67.60
D. F24880-8503	DUO Covered Station - 500ml (16oz) pH Neutralizing and 1000ml (32oz) 0.9% Saline, 12 x 14"	78.00




Sample Thief
 Hand Controlled
 Liquid Sampler

Polyethylene sampler has tapered tips welded to a 1 1/8" diameter tube, 36" long. End openings are 3/16" diameter. Overall length 41 1/2". Maximum capacity approximately 475ml. Other sizes made to order.

Catalog No.	List Price/Each
F37878-0000	\$93.00

Safety Acid/Caustic Samplers
 Retrieve Samples from Large Volume Tanks

- Use to sample sodium hydroxide, potassium hydroxide, hydrochloric acid, acetic acid, formaldehyde, sulfuric acid etc.
- Excellent for retrieval from tank cars, tank trucks and barges
- 500cc (16oz) polyethylene bottle screws onto a 4" long 1 1/2" diameter bar which has a slot for fluid entry
- 3ft. and 8ft. handle lengths available



Catalog No.	Description	List Price/Each
H37801-0000	(3ft.) handle (one-piece)	\$151.35
H37802-0000	(8ft.) handle (two-piece screw joint)	228.90

Battery Powered Engraving Tool
 Permanently Engraves in Seconds

Engrave test tubes, bottles, racks, keys, equipment, tools, electronics and more. Simply hold the Engraving Tool like a pen, press the power button, apply light pressure, and 'write' to permanently mark most materials including glass, plastic, wood, metal, bone, and stone.



Compact polypropylene stand (sold separately) holds the Battery Powered Engraving Tool and extra tip when not in use.

- Portable, cordless battery power; operates on two 'AA' batteries (not included)
- Engraving is permanent; marks will not fade if exposed to chemicals, sunlight or extreme temperatures
- Includes extra engraving tip

Catalog No.	Description	List Price/Each
F44150-0005	Battery Powered Engraving Tool, 7"L	\$21.00
F44150-0015	Engraving Tool Stand	11.50

Sterileware® Sampling Spatula
 V-shaped for easy retrieval of caked granules or powders



- Retrieve samples from bags and containers or directly from production lines
- Certified sterile, lot traceable, ready for one time use and disposal
- 3/4" x 6" spatula; 9" total length

Catalog No.	Description	List Price/Pack
H36930-0000	Sterile, polystyrene - 100 individually wrapped pieces	\$111.00/100
H37840-0000	Non-sterile, polypropylene - 12 pieces	21.80/12

Sterileware® Scoops
 Individually wrapped, sterilized and ready to go



- High-impact, white or red polystyrene scoops in popular sizes; high visibility red is ideal for light color samples
- Use to retrieve and weigh; flat bottom keeps them upright
- Certified sterile, lot traceable, ready for one time use and disposal
- Case of 100 or convenient 10 packs

Catalog No.	Color	Capacity	Bowl (L x W x H)	Total Length	List Price/Pack
H36902-0000	White	60ml (2oz)	3 x 1 x 1"	8"	\$105.75/100
H36902-0010	White	60ml (2oz)	3 x 1 x 1"	8"	20.85/10
H36904-0000	White	125ml (4oz)	3 1/2 x 1 x 1 3/4"	8 1/2"	\$135.75/100
H36904-0010	White	125ml (4oz)	3 1/2 x 1 x 1 3/4"	8 1/2"	26.25/10
H36906-0000	White	250ml (8oz)	5 x 2 x 2"	7 1/2"	\$181.65/100
H36906-0010	White	250ml (8oz)	5 x 2 x 2"	7 1/2"	35.10/10
H36902-2002	Red	60ml (2oz)	3 x 1 x 1"	8"	\$134.25/100
H36902-2010	Red	60ml (2oz)	3 x 1 x 1"	8"	26.25/10
H36904-2004	Red	125ml (4oz)	3 1/2 x 1 x 1 3/4"	8 1/2"	\$147.60/100
H36904-2010	Red	125ml (4oz)	3 1/2 x 1 x 1 3/4"	8 1/2"	28.35/10
H36906-2008	Red	250ml (8oz)	5 x 2 x 2"	7 1/2"	\$199.50/100
H36906-2010	Red	250ml (8oz)	5 x 2 x 2"	7 1/2"	39.45/10

