

## Soft-Grip Spatulas

The soft, silicone grip provides improved comfort and handling. The blade is constructed of stainless steel with a flexible, rounded end.

The handle is manufactured of polypropylene and layered with a silicone grip.

Handle dimensions are 4-1/8”L x 1”W x 1/2”H. Latex-free. Autoclavable.



Item No.	Overall Length	Blade Length	Blade Width	Quantity per Pack
SSSG03	7.25”	3”	3/4”	1
SSSG04	8.25”	4”	3/4”	1
SSSG05	9.25”	5”	3/4”	1
SSSG06	10.25”	6”	1”	1
SSSG08	12.25”	8”	1”	1

8221 East FM 917 • Alvarado, TX 76009  
 Phone: (800)369-5524 • Fax: (817)783-3571  
[sales@foxscientific.com](mailto:sales@foxscientific.com) • [www.foxscientific.com](http://www.foxscientific.com)