




### Economical, Reusable & Moisture Resistant Plastic

True north® freezer boxes by Heathrow Scientific, were designed with sustainability in mind and are made from a specially formulated corrugated or thin film polypropylene. True north freezer boxes are moisture resistant, which significantly reduces the chance of mold growth. With superior durability, polypropylene outlasts the cardboard box alternatives that are susceptible to softening and producing dust and can freeze (down to temperatures of -80°C) and thaw multiple times without a problem .

When it comes to saving valuable lab space, true north® Flatpack Freezer Boxes are supplied “flat-packed”, which eliminates the hassle of stacking and reduces the amount of space needed for storage. When ready to use, assembly is easy and pre-assembled dividers are included with each box. True north freezer boxes come complete in packs of 10, and multiple sizes and colors available.



**Freezer boxes are supplied “flat-packed”**  
Eliminating the hassle of stacking and reduces the amount of space needed for storage



Corrugated polypropylene 0.2/0.5/1.5-2.0/5 mL



Part No.	Tube	Wells	Well Dia.	I x w x h	Pkg.	Color
HS120379	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Natural <input type="checkbox"/>
HS120380	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120381	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120382	0.2 mL	144	7.9 mm	5.2 x 5.1 x 1.2 in (13.3 x 13.0 x 3.0 cm)	10/pk	Purple <input checked="" type="checkbox"/>
HS120383	0.5 mL	100	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Natural <input type="checkbox"/>
HS120384	0.5 mL	100	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120385	0.5 mL	100	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120386	0.5 mL	100	12 mm	5.2 x 5.1 x 1.6 in (13.3 x 13.0 x 4.2 cm)	10/pk	Purple <input checked="" type="checkbox"/>
HS120262	1.5/2.0 mL	81	12 mm	5.2 x 5.2 x 1.6 in (13.3 x 13.3 x 5.1 cm)	10/pk	Natural <input type="checkbox"/>
HS120376	1.5/2.0 mL	81	12 mm	5.2 x 5.2 x 1.6 in (13.3 x 13.3 x 5.1 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120377	1.5/2.0 mL	81	12 mm	5.2 x 5.2 x 1.6 in (13.3 x 13.3 x 5.1 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120378	1.5/2.0 mL	81	12 mm	5.2 x 5.2 x 1.6 in (13.3 x 13.3 x 5.1 cm)	10/pk	Purple <input checked="" type="checkbox"/>
HS120387	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.3 x 7.5 cm)	10/pk	Natural <input type="checkbox"/>
HS120388	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.3 x 7.5 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120389	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.3 x 7.5 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120390	5 mL	25	22.6 mm	5.2 x 5.1 x 2.9 in (13.3 x 13.3 x 7.5 cm)	10/pk	Purple <input checked="" type="checkbox"/>

Corrugated polypropylene 15 mL and 50 mL



Part No.	Tube	Wells	Well Dia.	I x w x h	Pkg.	Color
HS120391	15 mL	36	19.8 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Natural <input type="checkbox"/>
HS120392	15 mL	36	19.8 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120393	15 mL	36	19.8 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120394	15 mL	36	19.8 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Purple <input checked="" type="checkbox"/>
HS120366	50 mL	16	31.6 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Natural <input type="checkbox"/>
HS120367	50 mL	16	31.6 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120368	50 mL	16	31.6 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Red <input checked="" type="checkbox"/>
HS120369	50 mL	16	31.6 mm	5.7 x 5.78 x 4.72 in (14.5 x 14.7 x 12.0 cm)	10/pk	Purple <input checked="" type="checkbox"/>

Thin polypropylene film 1.5-2 mL



Part No.	Tube	Wells	Well Dia.	I x w x h	Pkg.	Color
HS120370	1.5/2.0 mL	81	12.1 mm	4.92 x 2.27 x 1.92 in (12.5 x 12.5 x 4.9 cm)	10/pk	White <input type="checkbox"/>
HS120371	1.5/2.0 mL	81	12.1 mm	4.92 x 2.27 x 1.92 in (12.5 x 12.5 x 4.9 cm)	10/pk	Blue <input checked="" type="checkbox"/>
HS120372	1.5/2.0 mL	81	12.1 mm	4.92 x 2.27 x 1.92 in (12.5 x 12.5 x 4.9 cm)	10/pk	Lt. Green <input checked="" type="checkbox"/>
HS120373	1.5/2.0 mL	81	12.1 mm	4.92 x 2.27 x 1.92 in (12.5 x 12.5 x 4.9 cm)	10/pk	Purple <input checked="" type="checkbox"/>



When finished with sample manipulation and ready for cold temperature storage, true north® freezer racks offer a variety of options to safely store samples in your freezer. You can choose from a wide selection of corrosion resistant stainless steel racks that are compatible with major brands or upright and chest freezers.

- All have handles for easy use
- Racks store boxes up to a 135 x 135 mm footprint and accommodate 2" tall boxes (up to 53 mm) as well as 3" tall boxes (up to 78 mm) depending on the rack mod



Fox Scientific, Inc. | 8221 East FM 917, Alvarado, TX 76009 USA  
 Tel.: +1-800-369-5524 | Fax: +1-817-783-3571 | E-mail: sales@foxscientific.com  
 Web: www.foxscientific.com