



Whatman®

934-AH™ RTU Glass Microfiber Filters

The 934-AH RTU contains the trusted binder-free, high purity borosilicate glass microfiber filter mentioned in Standard Methods. Now it comes pre-washed and pre-weighed. This eliminates filter pre-treating and streamlines your laboratory process. Each RTU filter is packaged within an aluminum pan with the filter weight clearly printed on a heat resistant label.

Features and Benefits

- ✦ **Fast and convenient:** Use the pre-washed, pre-weighed filter immediately – no need to pre-treat.
- ✦ **Efficient:** Superior fine particle retention at high flow rates.
- ✦ **High loading capacity:** Capable of processing very turbid samples.
- ✦ **High temperature tolerant:** Can be used with temperatures up to 550°C.



| Part Number | Qty | List Price | Promo Price |
|-------------|--------|------------|-------------|
| F100-001 | 100/pk | \$78.63 | \$66.84 |

Price good through April 31, 2012

ENVIRONMENTAL EXPRESS

ProWeigh Filters for TSS



All ProWeigh filters are prepared according to *Standard Methods* requirements. 47mm binderless glass fiber filters are weighed, washed 3 times in deionized water, dried at 105C and weighed to 0.1mg. Filters come packaged in a labeled aluminum pan with the weight of the filter printed on the label. In addition, the filter weight and ID number are printed in code 128 bar coding.

| Part Number | Qty | List Price |
|-------------|--------|------------|
| F895-001 | 100/pk | \$80.00 |



Fox Scientific, Inc.

8221 E. FM 917 • Alvarado, TX 76009
800-369-5524 Phone • 817-783-3571 Fax
www.foxscientific.com