



a xylem brand

3 REASONS TO UPGRADE YOUR BOD METER AND PROBE

EPA
Approved
Methodology

1

No more troublesome BOD blanks! The new MultiLab optical OBOD is accurate to +/- 0.1 mg/L. That's twice as accurate as a membrane probe making it easier than ever to meet your blank requirements.

2

Save time, money, and frustration with less maintenance. No more membranes with the MultiLab OBOD. The optical sensor cap needs replaced just once every 15 to 18 months.

3

Simplify your lab with one instrument. One MultiLab can measure three parameters at once - BOD, pH, Ammonia and more for compliance reporting.



Fox Scientific, Inc.

800-369-5524

sales@foxscientific.com

www.foxscientific.com

xylem
Let's Solve Water