


# Get Glass for your Clunker!

## Have a tired old hot plate/stirrer taking up space in your lab?

Retire your old worn out unit (regardless of make or model) and replace it with a new Corning® Hot Plate/Stirrer. Corning will reward you with free PYREX® Reusable Glass!



Select the model that most closely meets your needs!



### BENCH TOP HEATING

Hot plates available in 3 sizes including digital units with LED temperature displays.

	Top Plate Size
<b>Digital:</b>	
PC 400D (6795-400D)	5" x 7"
PC 600D (6795-600D)	10" x 10"
<b>Analog:</b>	
PC200 (6795-200)	4" x 5"



### BENCH TOP STIRRING

Single position stirrers available in 4 sizes including digital units with LED stirring speed displays.

	Top Plate Size
<b>Digital:</b>	
PC 410D (6795-410D)	5" x 7"
PC 610D (6795-610D)	10" x 10"
<b>Analog:</b>	
PC 210 (6795-210)	4" x 5"
PC 611 (6795-611)	11" x 11"



### COMBINED STIRRING AND HEATING

Single position stirring hot plates available in 3 sizes including digital units with LED stirring speed and temperature displays.

	Top Plate Size
<b>Digital:</b>	
PC 420D (6795-420D)	5" x 7"
PC 620D (6795-620D)	10" x 10"
<b>Analog:</b>	
PC 220 (6795-220)	4" x 5"

For detailed information on Corning's hot plates and stirrers, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**CORNING**

### Corning® Digital Hot Plates/Stirrers Offer:

- ▶ Digital LED temperature display is adjustable in 5°C increments. Display blinks until set temperature is reached.
- ▶ Microprocessor maintains consistent and repeatable temperature settings from 5°C (if ambient temperature is 0°C or lower) up to 550°C.

### Digital Hot Plate Safety Features:

- ▶ Bright LED hot top icon lights up when top plate is over 60°C, even when heat control is turned off.
- ▶ Separate temperature sensor provides power cutoff if unit overheats.

### Get Free PYREX® Glassware! Participation is Easy:

- ▶ Purchase a designated Corning hot plate, stirrer, or combination unit between November 1, 2009 and January 31, 2010. Offer good for up to five units.
- ▶ Complete the information requested below and return, with your proof of purchase (invoice or packing slip) by one of the following methods:
  1. Email: [lifesciences@corning.com](mailto:lifesciences@corning.com)
  2. Fax: **978.442.2478**
  3. Mail to: Attn: Glass for Clunker Offer, Corning Incorporated, Life Sciences, Tower 2, 4th Floor, 900 Chelmsford St., Lowell, MA 01851

### Purchase one or more hot plate/stirrer models and apply for the appropriate GLASS reward:

	SELECT
▶ Models PC 200, 210, 400D, 410D, receive <b>2</b> cases of glassware	<input type="checkbox"/>
▶ Models PC 220, 600D, 420D, 610D, receive <b>3</b> cases of glassware	<input type="checkbox"/>
▶ Models PC 611, 620D, receive <b>4</b> cases of glassware	<input type="checkbox"/>

### Select Your Free Glass Product:



**BEAKERS:** Cat. No's. 1000-50, 100, 150, 250, 400, 600, 1L

**STORAGE BOTTLES:** Cat. No's. 1395-100, 250, 500, 1L, 5L

**VOLUMETRIC FLASKS:** Cat. No's. 5640-5, 10, 25, 50, 100, 250, 500, 1L, 2L

**ERLENMEYER FLASK WITH SCREW CAP:** Cat. No's. 4995-250, 500, 1L, 2L, 4L, 6L

Indicate number of cases and catalog number of free product to be shipped:

Unit you will be retiring:

Brand \_\_\_\_\_ Model \_\_\_\_\_ Age \_\_\_\_\_

Your Name \_\_\_\_\_

Title \_\_\_\_\_

Company \_\_\_\_\_

Address (No P.O. box, please.) \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip + 4 \_\_\_\_\_

Telephone \_\_\_\_\_ Fax \_\_\_\_\_

E-mail \_\_\_\_\_

Unit does not need to be returned to Corning.

Offer valid on purchases made between November 1, 2009 and January 31, 2010. All requests must be postmarked no later than February 26, 2010. Offer void where prohibited by government or company policy, or local law. Offer valid in the U.S. and Canada only.

# CORNING

Corning Incorporated  
Life Sciences

Tower 2, 4th Floor  
900 Chelmsford St.  
Lowell, MA 01851  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)