

Thermo Scientific Orion Star A221

Advanced pH portable meter

The Thermo Scientific™ Orion™ Star A221 pH Portable Meter makes necessary testing easier. Its durability and dependability make it a great choice for samples every day, anywhere.



Get the readings that you need where you need them with the Thermo Scientific Orion Star A221 pH Portable Meter. The large, backlit graphic LCD display shows pH/mV/ORP and temperature along with electrode status, time, date and calibration points.

Features and Benefits

- Easy meter operation assured with on-screen text prompts, menu-specific function keys and a multi-language interface for English, Spanish, French, Italian, German and Chinese
- Don't miss a reading—AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when continuous readings are stable
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers or manual option for custom buffer values
- Non-volatile memory holds up to 1000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary Orion Star Com computer software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately)
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty

Fox
Scientific, Inc.
Thermo
SCIENTIFIC

Thermo Scientific Orion Star A221 pH Portable Meter

Specifications	
pH	Range: -2.000 to 20.000 Resolution: 0.1, 0.01, 0.001 Relative Accuracy: ± 0.002 Calibration Points: Up to 5
mV/RmV	Range: ± 2000.0 mV Resolution: 0.1 mV Relative Accuracy: ± 0.2 mV or ± 0.05 % of reading whichever is greater E_H ORP Mode: Yes
Temperature	Range: -5 to 105 °C, 23 to 221 °F Resolution: 0.1 °C, 0.1 °F Relative Accuracy: ± 0.1 °C, ± 0.1 °F Offset Calibration: 1 point
Datalogging	Memory: 1000 with time and date stamp Log Function: Manual, Automatic with AUTO-READ
Inputs	pH Electrode: BNC, reference pin ATC Probe: 8-pin mini-DIN
Output	RS232, USB
Power	AC Adapter: Optional—universal, 100-240 VAC Battery Power: Included—4 AAs Battery Life: 800 hrs

Meter Ordering Information	
Part No.	Description
STARA2210	Orion Star A221 pH Portable Meter - Four AA batteries (installed), literature, PC cable and test certificate
STARA2215	Orion Star A221 pH Portable Meter Standard Kit - 8107UWMMMD Orion ROSS Ultra Triode gel epoxy pH/ATC electrode, 3 m cable - Orion pH 4, 7, 10 buffer and rinse pouches, 10 each - Orion ROSS pH electrode storage solution, 475 mL - Protective meter armor with electrode holder - Hard-sided field case - Four AA batteries (installed), literature, PC cable and test certificate

STARA2216

Orion Star A221 pH Portable Meter

General Purpose Kit

- 9107BNMD Orion Triode gel epoxy pH/ATC electrode
- Orion pH 4, 7, 10 buffer and rinse pouches, 10 each
- Orion pH electrode storage solution, 475 mL
- Protective meter armor with electrode holder
- Hard-sided field case
- Four AA batteries (installed), literature, PC cable and test certificate

STARA2217

Orion Star A221 pH Portable Meter ORP Kit

- 9179BNMD Orion gel epoxy ORP/ATC electrode
- 967961 Orion ORP standard, 5 x 60 mL
- 911110 Orion rinse solution pouch 10 pack
- 910060 Orion storage solution, 5 x 60 mL
- Protective meter armor with electrode holder
- Hard-sided field case
- Four AA batteries (installed), literature, PC cable and test certificate

STARA2218

Orion Star A221 pH Portable Meter

Food Sample Kit

- 8163BNWP Orion ROSS Spear-Tip pH electrode
- 927007MD Orion stainless steel ATC probe
- Orion pH 4, 7, 10 buffer and rinse pouches, 10 each
- ROSS pH electrode storage solution, 475 mL
- Protective meter armor with electrode holder
- Hard-sided field case
- Four AA batteries (installed), literature, PC cable and test certificate

Accessories Ordering Information

Part No.	Description
STARA-CS	Orion Star A series portable meter hard-sided field case
1210005	Orion Star A series portable meter soft-sided field case
STARA-AR	Orion Star A series portable meter armor with pH electrode holder, conductivity and DO probe holder, stand and adjustable hand-strap
STARA-ESPH	Orion pH electrode holder for portable Star A series meter armor
810017	Orion pH electrode storage sleeve for pH electrode holder
1010003	Universal power adapter, 100-240 V, 50/60 Hz

For more information, contact your local Thermo Scientific products dealer or call our customer and technical service experts at 1-800-225-1480 (for the US and Canada) or visit www.thermoscientific.com/water.



Fox Scientific, Inc.

8221 E FM 917 Alvarado, TX 76009

800-369-5524

sales@foxscientific.com

www.foxscientific.com

thermoscientific.com/water

© 2016 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Water and Lab Products

North America

Toll Free: 1-800-225-1480

Tel: 1-978-232-6000

info.water@thermofisher.com

Germany

Tel: (49) 6184-90-6000

info.water.uk@thermofisher.com

China

Tel: (86) 21-68654588

wai.asia@thermofisher.com

India

Tel: (91) 22-4157-8800

wai.asia@thermofisher.com

Singapore

Tel: (65) 6778-6876

wai.asia@thermofisher.com

Japan

Tel: (81) 045-453-9175

wai.asia@thermofisher.com

Australia

Tel: (613) 9757-4300

In Australia: (1300) 735-295

InfoWaterAU@thermofisher.com

S-STARA221-E 0216 RevD



Fox
Scientific, Inc.
Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand