

*Fox Scientific, Inc.*



# Clearance Pricing

**4600-0121C 120-mL sterile coliform bottle, with thiosulfate, screw cap, tamper seal, (polystyrene), 100/pk**

- Non-fluorescent Polystyrene
- 120mL Bottle with 100mL fill line
- Fill Line accuracy of 100mL +/- 2.5%
- Screw cap with tamper-proof seal
- Contains 50mg sodium thiosulfate tablet certified to neutralize a minimum of 10ppm chlorine
- Lot number and expiration date laser marked on each bottle
- Lot specific Certificate of Conformance included with each shipment

Our sterile coliform sample bottles have a 120-mL volume to meet EPA requirements, with a 100-mL mark molded into the polystyrene plastic, encircling the bottle to make it easy to see and feel the fill line. The bottle contains a sodium thiosulfate tablet to neutralize chlorine and has a tamper proof seal to ensure sterilization. Sold as a box of 100 bottles.

(800) 369-5524  
sales@foxscientific.com  
www.foxscientific.com



**\$87.51/PK**

Price valid for current inventory only!